

For Research Use Only

NDUFS6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68329-6-PBS



Basic Information

Catalog Number:

68329-6-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6286

GenBank Accession Number:

BC046155

GeneID (NCBI):

4726

UNIPROT ID:

O75380

Full Name:

NADH dehydrogenase (ubiquinone)

Fe-S protein 6, 13kDa (NADH-coenzyme Q reductase)

Calculated MW:

14 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C11A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50662-3, NDUFS6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68329-4-PBS. Detection antibody: 68329-6-PBS. Standard: Ag6286. Range: 0.098-1.563 ng/mL.