

For Research Use Only

GPKOW Monoclonal antibody, PBS Only



Catalog Number: 68338-1-PBS

Basic Information

Catalog Number: 68338-1-PBS	GenBank Accession Number: BC000397	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 27238	CloneNo.: 1G11C4
Source: Mouse	UNIPROT ID: Q92917	
Isotype: IgG1	Full Name: G patch domain and KOW motifs	
Immunogen Catalog Number: AG6545	Calculated MW: 52 kDa	
	Observed MW: 52 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat

Background Information

GPKOW is a RNA-binding protein involved in pre-mRNA splicing. It interacts directly with protein kinase A and protein kinase X and is also found associated with the spliceosome.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68338-1-Ig (GPKOW antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68338-1-PBS in a different storage buffer formulation.