

For Research Use Only

NDUFB7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68362-3-PBS



Basic Information

Catalog Number:

68362-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG6720

GenBank Accession Number:

BC002595

GeneID (NCBI):

4713

UNIPROT ID:

P17568

Full Name:

NADH dehydrogenase (ubiquinone) 1
beta subcomplex, 7, 18kDa

Calculated MW:

16 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2B1C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50699-1, NDUFB7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68362-2-PBS. Detection antibody: 68362-3-PBS. Standard: Ag6720. Range: 0.098-100 ng/mL.