

For Research Use Only

DSP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68364-3-PBS



Basic Information

Catalog Number: 68364-3-PBS	GenBank Accession Number: BC140802	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1832	CloneNo.: 2D5C2
Source: Mouse	UNIPROT ID: P15924	
Isotype: IgG1	Full Name: desmoplakin	
Immunogen Catalog Number: AG30499	Calculated MW: 2871 aa, 332 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51116-1, DSP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68364-2-PBS. Detection antibody: 68364-3-PBS. Standard: Ag30499. Range: 0.195-50 ng/mL