

For Research Use Only

MITD1 Monoclonal antibody, PBS Only



Catalog Number: 68367-1-PBS

Basic Information

Catalog Number: 68367-1-PBS	GenBank Accession Number: BC018453	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 129531	CloneNo.: 1F11E2
Source: Mouse	UNIPROT ID: Q8WV92	
Isotype: IgG1	Full Name: MIT, microtubule interacting and transport, domain containing 1	
Immunogen Catalog Number: AG11237	Calculated MW: 249 aa, 29 kDa	
	Observed MW: 30 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, Mouse, Rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68367-1-Ig (MITD1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Beta Actin Monoclonal antibody (HRP-66009) as loading control. This data was developed using the same antibody clone with 68367-1-PBS in a different storage buffer formulation.