

For Research Use Only

CREB3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68379-6-PBS



Basic Information

Catalog Number: 68379-6-PBS	GenBank Accession Number: BC009402	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10488	CloneNo.: 3E1C10
Source: Mouse	UNIPROT ID: O43889	
Isotype: IgG1	Full Name: cAMP responsive element binding protein 3	
Immunogen Catalog Number: AG1823	Calculated MW: 44 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50828-3, CREB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68379-5-PBS. Detection antibody: 68379-6-PBS. Standard: Ag1823. Range: 0.195-100 ng/mL.