

For Research Use Only

RRP8 Monoclonal antibody

Catalog Number: 68382-1-Ig



Basic Information

Catalog Number:

68382-1-Ig

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16803

GenBank Accession Number:

BC001071

GeneID (NCBI):

23378

UNIPROT ID:

O43159

Full Name:

ribosomal RNA processing 8,
methyltransferase, homolog (yeast)

Calculated MW:

456 aa, 51 kDa

Observed MW:

51 kDa

Purification Method:

Protein G purification

CloneNo.:

1F5C2

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse, rat

Positive Controls:

WB : U2OS cells, LNCaP cells, HeLa cells, HEK-293 cells,
Jurkat cells, K-562 cells, HSC-T6 cells, 4T1 cells,
NIH/3T3 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68382-1-Ig (RRP8 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control.