

For Research Use Only

RRP8 Monoclonal antibody, PBS Only



Catalog Number: 68382-1-PBS

Basic Information

Catalog Number:

68382-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16803

GenBank Accession Number:

BC001071

GeneID (NCBI):

23378

UNIPROT ID:

O43159

Full Name:

ribosomal RNA processing 8,
methyltransferase, homolog (yeast)

Calculated MW:

456 aa, 51 kDa

Observed MW:

51 kDa

Purification Method:

Protein G purification

CloneNo.:

1F5C2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68382-1-Ig (RRP8 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68382-1-PBS in a different storage buffer formulation.