

For Research Use Only

TNIP1 Monoclonal antibody, PBS Only



Catalog Number: 68387-1-PBS

Basic Information

Catalog Number: 68387-1-PBS	GenBank Accession Number: BC014008	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10318	CloneNo.: 3F1E6
Source: Mouse	UNIPROT ID: Q15025	
Isotype: IgG1	Full Name: TNFAIP3 interacting protein 1	
Immunogen Catalog Number: AG8885	Calculated MW: 72 kDa	
	Observed MW: 72 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Background Information

TNIP1, also named as KIAA0113, NAF1, Nip40-1, VAN and ABIN-1, is involved in inhibition of nuclear factor- B (NF- B) activation by interacting with TNF alpha-induced protein 3, A20 (TNFAIP3). TNIP1 repression of RARs in the presence of RA places it in a small category of corepressors of agonist-bound NRs.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS only

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68387-1-Ig (TNIP1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68387-1-PBS in a different storage buffer formulation.