

For Research Use Only

RFC4 Monoclonal antibody, PBS Only

Catalog Number: 68390-1-PBS



Basic Information

Catalog Number:

68390-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33295

GenBank Accession Number:

BC017452

GeneID (NCBI):

5984

UNIPROT ID:

P35249

Full Name:

replication factor C (activator 1) 4,

37kDa

Calculated MW:

37 kDa

Observed MW:

37-40 kDa

Purification Method:

Protein G purification

CloneNo.:

1A6A11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68390-1-Ig (RFC4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 68390-1-PBS in a different storage buffer formulation.