

For Research Use Only

ERGIC3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68398-2-PBS



Basic Information

Catalog Number: 68398-2-PBS	GenBank Accession Number: BC009765	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 51614	CloneNo.: 1D12B11
Source: Mouse	UNIPROT ID: Q9Y282	
Isotype: IgG1	Full Name: ERGIC and golgi 3	
Immunogen Catalog Number: AG9190	Calculated MW: 383 aa, 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50997-1, ERGIC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68398-2-PBS. Detection antibody: 68398-3-PBS. Standard:Ag9190. Range: 0.195-100 ng/mL.