

For Research Use Only

# ERGIC3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68398-4-PBS



## Basic Information

Catalog Number: 68398-4-PBS	GenBank Accession Number: BC009765	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51614	CloneNo.: 1D8H5
Source: Mouse	UNIPROT ID: Q9Y282	
Isotype: IgG2a	Full Name: ERGIC and golgi 3	
Immunogen Catalog Number: AG9190	Calculated MW: 383 aa, 43 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50997-2, ERGIC3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68398-4-PBS. Detection antibody: 68398-3-PBS. Standard:Ag9190. Range: 0.195-100 ng/mL.