

For Research Use Only

GATC Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68405-2-PBS



Basic Information

Catalog Number:

68405-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG19857

GenBank Accession Number:

BC107145

GeneID (NCBI):

283459

UNIPROT ID:

O43716

Full Name:

glutamyl-tRNA(Gln)

amidotransferase, subunit C homolog
(bacterial)

Calculated MW:

136 aa, 15 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1H6F11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50658-1, GATC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68405-2-PBS. Detection antibody: 68405-3-PBS. Standard: Ag19857. Range: 0.098-3.125 ng/mL.