

For Research Use Only

# GATC Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68405-4-PBS



## Basic Information

Catalog Number:

68405-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19857

GenBank Accession Number:

BC107145

GeneID (NCBI):

283459

UNIPROT ID:

O43716

Full Name:

glutamyl-tRNA(Gln)  
amidotransferase, subunit C homolog  
(bacterial)

Calculated MW:

136 aa, 15 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C2H12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50658-2, GATC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68405-4-PBS. Detection antibody: 68405-5-PBS. Standard: Ag19857. Range: 0.781-100 ng/mL.