

For Research Use Only

GATC Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68405-5-PBS



Basic Information

Catalog Number: 68405-5-PBS	GenBank Accession Number: BC107145	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 283459	CloneNo.: 1D1H11
Source: Mouse	UNIPROT ID: O43716	
Isotype: IgG2a	Full Name: glutamyl-tRNA(Gln) amidotransferase, subunit C homolog (bacterial)	
Immunogen Catalog Number: AG19857	Calculated MW: 136 aa, 15 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50658-2, GATC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68405-4-PBS. Detection antibody: 68405-5-PBS. Standard: Ag19857. Range: 0.781-100 ng/mL.