

For Research Use Only

# TBC1D23 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68409-2-PBS



## Basic Information

Catalog Number: 68409-2-PBS	GenBank Accession Number: BC010221	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 55773	CloneNo.: 2C11G11
Source: Mouse	UNIPROT ID: Q9NUY8	
Isotype: IgG1	Full Name: TBC1 domain family, member 23	
Immunogen Catalog Number: AG10727	Calculated MW: 699 aa, 78 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51051-1, TBC1D23 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68409-2-PBS. Detection antibody: 68409-3-PBS. Standard:Ag13334. Range: 6.25-200 ng/mL