

For Research Use Only

HSDL2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68414-5-PBS



Basic Information

| | | |
|-------------------------------------|---|--|
| Catalog Number: 68414-5-PBS | GenBank Accession Number: BC004331 | Purification Method: Protein G purification |
| Size: 1 mg/ml | GeneID (NCBI): 84263 | CloneNo.: 3A9C11 |
| Source: Mouse | UNIPROT ID: Q6YN16 | |
| Isotype: IgG1 | Full Name: hydroxysteroid dehydrogenase like 2 | |
| Immunogen Catalog Number: AG8314 | Calculated MW: 45 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50441-3, HSDL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68414-5-PBS. Detection antibody: 68414-6-PBS. Standard: Ag8314. Range: 0.098-100 ng/mL.