

For Research Use Only

MRPL38 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68418-3-PBS



Basic Information

Catalog Number: 68418-3-PBS	GenBank Accession Number: BC013311	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 64978	CloneNo.: 1C3G2
Source: Mouse	UNIPROT ID: Q96DV4	
Isotype: IgG1	Full Name: mitochondrial ribosomal protein L38	
Immunogen Catalog Number: AG8927	Calculated MW: 346 aa, 41 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50538-1, MRPL38 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68418-2-PBS. Detection antibody: 68418-3-PBS. Standard: Ag8927. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50538-4, MRPL38 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68418-1-PBS. Detection antibody: 68418-3-PBS. Standard: Ag8927. Range: 6.25-100 ng/mL.