

For Research Use Only

UBOX5 Monoclonal antibody, PBS Only



Catalog Number: 68422-1-PBS

Basic Information

Catalog Number: 68422-1-PBS	GenBank Accession Number: BC000515	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 22888	CloneNo.: 2A7A4
Source: Mouse	UNIPROT ID: O94941	
Isotype: IgG1	Full Name: U-box domain containing 5	
Immunogen Catalog Number: AG6573	Calculated MW: 59 kDa	
	Observed MW: 60 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68422-1-Ig (UBOX5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68422-1-PBS in a different storage buffer formulation.



Various lysates were subjected to SDS PAGE followed by western blot with 68422-1-Ig (UBOX5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68422-1-PBS in a different storage buffer formulation.