

For Research Use Only

# ZNF174 Monoclonal antibody, PBS Only



Catalog Number: 68426-1-PBS

## Basic Information

Catalog Number: 68426-1-PBS	GenBank Accession Number: BC000876	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7727	CloneNo.: 2D9D3
Source: Mouse	UNIPROT ID: Q15697	
Isotype: IgG3	Full Name: zinc finger protein 174	
Immunogen Catalog Number: AG12967	Calculated MW: 46 kDa	
	Observed MW: 46 kDa	

## Applications

Tested Applications:  
WB, Indirect ELISA

Species Specificity:  
Human, mouse, rat, rabbit

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68426-1-Ig (ZNF174 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68426-1-PBS in a different storage buffer formulation.