

For Research Use Only

# STX4 Monoclonal antibody, PBS Only

Catalog Number: 68429-1-PBS



## Basic Information

Catalog Number:

68429-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7053

GenBank Accession Number:

BC002436

GeneID (NCBI):

6810

UNIPROT ID:

Q12846

Full Name:

syntaxin 4

Calculated MW:

34 kDa

Observed MW:

35 kDa

Purification Method:

Protein G purification

CloneNo.:

4D6C3

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68429-1-Ig (STX4 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68429-1-PBS in a different storage buffer formulation.