

For Research Use Only

Synaptotagmin-9 Monoclonal antibody, PBS Only



Catalog Number: 68430-1-PBS

Basic Information

| | | |
|---|--|---|
| Catalog Number: 68430-1-PBS | GenBank Accession Number: BC029605 | Purification Method: Protein G purification |
| Size: 1 mg/ml | GeneID (NCBI): 143425 | CloneNo.: 3D2F2 |
| Source: Mouse | UNIPROT ID: Q86556 | |
| Isotype: IgG1 | Full Name: synaptotagmin IX | |
| Immunogen Catalog Number: AG33194 | Calculated MW: 491 aa, 52 kDa | |
| | Observed MW: 52 kDa | |

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat, rabbit

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68430-1-Ig (Synaptotagmin-9 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68430-1-PBS in a different storage buffer formulation.