

For Research Use Only

BHLHE40 Monoclonal antibody, PBS Only



Catalog Number: 68435-1-PBS

Basic Information

Catalog Number:

68435-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG12236

GenBank Accession Number:

BC082238

GeneID (NCBI):

8553

UNIPROT ID:

O14503

Full Name:

basic helix-loop-helix family,
member e40

Calculated MW:

412 aa, 46 kDa

Observed MW:

46-50 kDa

Purification Method:

Protein A purification

CloneNo.:

1A3E1

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

Background Information

BHLHE40 (Basic Helix-Loop-Helix Family Member E40), also known as BHLHB2, STRA13, DEC1, or SHARP2, is a member of the basic helix-loop-helix (bHLH) protein family, a large superfamily of transcriptional regulators expressed in many organisms. BHLHE40 is known to regulate a wide variety of essential cellular processes, including cell cycle, cellular proliferation, programmed cell death, cellular development and differentiation, as well as circadian rhythms (PMID: 3451158). It is reported that BHLHE40 is overexpressed in gastric, breast, and brain tumors; and downregulated in colorectal, esophageal, pancreatic and lung cancer (PMID: 32577154).

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68435-1-Ig (BHLHE40 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68435-1-PBS in a different storage buffer formulation.