

For Research Use Only

# FBP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68446-3-PBS



## Basic Information

Catalog Number: 68446-3-PBS	GenBank Accession Number: BC012927	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 2203	CloneNo.: 3G7C8
Source: Mouse	UNIPROT ID: P09467	
Isotype: IgG1	Full Name: fructose-1,6-bisphosphatase 1	
Immunogen Catalog Number: AG3837	Calculated MW: 338 aa, 37 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50428-1, FBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68446-2-PBS. Detection antibody: 68446-3-PBS. Standard: Ag3837. Range: 1.563-100 ng/mL.