

For Research Use Only

FBP1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68446-5-PBS



Basic Information

Catalog Number: 68446-5-PBS	GenBank Accession Number: BC012927	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 2203	CloneNo.: 1D8E5
Source: Mouse	UNIPROT ID: P09467	
Isotype: IgG1	Full Name: fructose-1,6-bisphosphatase 1	
Immunogen Catalog Number: AG3837	Calculated MW: 338 aa, 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50428-2, FBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68446-4-PBS. Detection antibody: 68446-5-PBS. Standard: Ag3837. Range: 3.125-100 ng/mL.