

For Research Use Only

OTUD6B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68450-4-PBS



Basic Information

Catalog Number:

68450-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33607

GenBank Accession Number:

BC029760

GeneID (NCBI):

51633

UNIPROT ID:

Q8N6M0

Full Name:

OTU domain containing 6B

Calculated MW:

293 aa, 34 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F11C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50962-3, OTUD6B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68450-1-PBS. Detection antibody: 68450-4-PBS. Standard: Ag33607. Range: 0.098-12.5 ng/mL