

For Research Use Only

GIMAP7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68451-4-PBS



Basic Information

Catalog Number: 68451-4-PBS	GenBank Accession Number: BC027613	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 168537	CloneNo.: 1C3C3
Source: Mouse	UNIPROT ID: Q8NHV1	
Isotype: IgG1	Full Name: GTPase, IMAP family member 7	
Immunogen Catalog Number: AG10740	Calculated MW: 300 aa, 35 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50958-3, GIMAP7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68451-1-PBS. Detection antibody: 68451-4-PBS. Standard:Ag10740. Range: 0.098-100 ng/mL.