

For Research Use Only

FBLN1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68453-2-PBS



Basic Information

Catalog Number: 68453-2-PBS	GenBank Accession Number: BC022497	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 2192	CloneNo.: 1G11E2
Source: Mouse	UNIPROT ID: P23142	
Isotype: IgG1	Full Name: fibulin 1	
Immunogen Catalog Number: AG14249	Calculated MW: 703 aa, 77 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50734-1, FBLN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68453-2-PBS. Detection antibody: 68453-3-PBS. Standard:Ag14249. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50734-3, FBLN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68453-2-PBS. Detection antibody: 68453-1-PBS. Standard:Ag14249. Range: 0.391-50 ng/mL.