

For Research Use Only

DDX43 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68454-2-PBS

Basic Information

Catalog Number:

68454-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11484

GenBank Accession Number:

BC066938

GeneID (NCBI):

55510

UNIPROT ID:

Q9NXZ2

Full Name:

DEAD (Asp-Glu-Ala-Asp) box polypeptide 43

Calculated MW:

648 aa, 73 kDa

Purification Method:

Protein G purification

CloneNo.:

1C3F9

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50126-1, DDX43 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68454-2-PBS. Detection antibody: 68454-3-PBS. Standard: Ag11484. Range: 6.25-100 ng/mL