

For Research Use Only

DDX43 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68454-3-PBS

Basic Information

Catalog Number: 68454-3-PBS	GenBank Accession Number: BC066938	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 55510	CloneNo.: 2A1H3
Source: Mouse	UNIPROT ID: Q9NXZ2	
Isotype: IgG1	Full Name: DEAD (Asp-Glu-Ala-Asp) box polypeptide 43	
Immunogen Catalog Number: AG11484	Calculated MW: 648 aa, 73 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50126-1, DDX43 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68454-2-PBS. Detection antibody: 68454-3-PBS. Standard: Ag11484. Range: 6.25-100 ng/mL