

For Research Use Only

EIF2S2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68463-2-PBS



Basic Information

Catalog Number:

68463-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18349

GenBank Accession Number:

BC000461

GeneID (NCBI):

8894

UNIPROT ID:

P20042

Full Name:

eukaryotic translation initiation
factor 2, subunit 2 beta, 38kDa

Calculated MW:

38 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B2D4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50725-1, EIF2S2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68463-2-PBS. Detection antibody: 68463-3-PBS. Standard: Ag18349. Range: 0.391-100 ng/mL.