

For Research Use Only

PPA2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68468-5-PBS



Basic Information

Catalog Number:

68468-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG10125

GenBank Accession Number:

BC057219

GeneID (NCBI):

27068

UNIPROT ID:

Q9H2U2

Full Name:

pyrophosphatase (inorganic) 2

Calculated MW:

305aa,35 kDa; 334aa,38 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F8G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50620-2, PPA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68468-4-PBS. Detection antibody: 68468-5-PBS. Standard: Ag10125. Range: 0.391-25 ng/mL