

For Research Use Only

IGF2BP3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68470-4-PBS



Basic Information

Catalog Number: 68470-4-PBS	GenBank Accession Number: BC065269	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10643	CloneNo.: 1B11D8
Source: Mouse	UNIPROT ID: O00425	
Isotype: IgG2a	Full Name: insulin-like growth factor 2 mRNA binding protein 3	
Immunogen Catalog Number: AG6349	Calculated MW: 64 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51098-2, IGF2BP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68470-4-PBS. Detection antibody: 68470-3-PBS. Standard:Ag6349. Range: 0.781-50 ng/mL