

For Research Use Only

# IGF2BP3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68470-5-PBS



## Basic Information

<b>Catalog Number:</b> 68470-5-PBS	<b>GenBank Accession Number:</b> BC065269	<b>Purification Method:</b> Protein A Magarose purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 10643	<b>CloneNo.:</b> 2F4F8
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> O00425	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> insulin-like growth factor 2 mRNA binding protein 3	
<b>Immunogen Catalog Number:</b> AG6349	<b>Calculated MW:</b> 64 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51098-3, IGF2BP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68470-5-PBS. Detection antibody: 68470-3-PBS. Standard: Ag6349. Range: 3.125-100 ng/mL.