

For Research Use Only

# MAOA Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68471-2-PBS



## Basic Information

|                                      |                                       |                                                         |
|--------------------------------------|---------------------------------------|---------------------------------------------------------|
| Catalog Number:<br>68471-2-PBS       | GenBank Accession Number:<br>BC008064 | Purification Method:<br>Protein G Magarose purification |
| Size:<br>1 mg/ml                     | GeneID (NCBI):<br>4128                | CloneNo.:<br>2D7E5                                      |
| Source:<br>Mouse                     | UNIPROT ID:<br>P21397                 |                                                         |
| Isotype:<br>IgG1                     | Full Name:<br>monoamine oxidase A     |                                                         |
| Immunogen Catalog Number:<br>AG30028 | Calculated MW:<br>60 kDa              |                                                         |

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50632-1, MAOA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68471-2-PBS. Detection antibody: 68471-3-PBS. Standard: Ag30028. Range: 0.098-100 ng/mL.