

For Research Use Only

MEN1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68473-3-PBS



Basic Information

Catalog Number: 68473-3-PBS	GenBank Accession Number: BC002544	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4221	CloneNo.: 3E4G10
Source: Mouse	UNIPROT ID: O00255	
Isotype: IgG1	Full Name: multiple endocrine neoplasia I	
Immunogen Catalog Number: AG7566	Calculated MW: 68 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50886-2, MEN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68473-1-PBS. Detection antibody: 68473-3-PBS. Standard: Ag7566. Range: 3.125-200 ng/mL.