

For Research Use Only

EBP2 Monoclonal antibody, PBS Only

Catalog Number: 68477-1-PBS



Basic Information

Catalog Number: 68477-1-PBS	GenBank Accession Number: BC009175	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10969	CloneNo.: 1F6G1
Source: Mouse	UNIPROT ID: Q99848	
Isotype: IgG1	Full Name: EBNA1 binding protein 2	
Immunogen Catalog Number: AG8833	Calculated MW: 306 aa, 35 kDa	
	Observed MW: 38-40 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human, mouse, rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68477-1-Ig (EBP2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68477-1-PBS in a different storage buffer formulation.