

For Research Use Only

EIF2S1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68479-2-PBS

Basic Information

Catalog Number: 68479-2-PBS	GenBank Accession Number: BC002513	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1965	CloneNo.: 2A11C1
Source: Mouse	UNIPROT ID: P05198	
Isotype: IgG1	Full Name: eukaryotic translation initiation factor 2, subunit 1 alpha, 35kDa	
Immunogen Catalog Number: AG28309	Calculated MW: 36 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50179-1, EIF2S1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68479-2-PBS. Detection antibody: 68479-3-PBS. Standard: Ag28309. Range: 0.195-25 ng/mL