

For Research Use Only

# EIF2S1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68479-3-PBS

## Basic Information

Catalog Number:

68479-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28309

GenBank Accession Number:

BC002513

GeneID (NCBI):

1965

UNIPROT ID:

P05198

Full Name:

eukaryotic translation initiation factor 2, subunit 1 alpha, 35kDa

Calculated MW:

36 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B6C8

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50179-1, EIF2S1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68479-2-PBS. Detection antibody: 68479-3-PBS. Standard: Ag28309. Range: 0.195-25 ng/mL