

For Research Use Only

DIABLO Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68480-2-PBS

Basic Information

Catalog Number: 68480-2-PBS	GenBank Accession Number: BC011909	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 56616	CloneNo.: 1H4C8
Source: Mouse	UNIPROT ID: Q9NR28	
Isotype: IgG2a	Full Name: diablo homolog (Drosophila)	
Immunogen Catalog Number: AG32774	Calculated MW: 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50334-1, DIABLO Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68480-2-PBS. Detection antibody: 68480-3-PBS. Standard: Ag32774. Range: 0.098-3.125 ng/mL.