

For Research Use Only

# DIABLO Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68480-5-PBS

## Basic Information

Catalog Number:

68480-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG32774

GenBank Accession Number:

BC011909

GeneID (NCBI):

56616

UNIPROT ID:

Q9NR28

Full Name:

diablo homolog (Drosophila)

Calculated MW:

27 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A2A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50334-3, DIABLO Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68480-4-PBS. Detection antibody: 68480-5-PBS. Standard: Ag32774. Range: 0.781-50 ng/mL