

For Research Use Only

EIF2AK3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68482-2-PBS



Basic Information

Catalog Number:

68482-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG18254

GenBank Accession Number:

BC126354

GeneID (NCBI):

9451

UNIPROT ID:

Q9NZJ5

Full Name:

eukaryotic translation initiation

factor 2-alpha kinase 3

Calculated MW:

1116 aa, 125 kDa

Purification Method:

Protein A purification

CloneNo.:

4B3E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50733-1, EIF2AK3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68482-2-PBS. Detection antibody: 68482-3-PBS. Standard:Ag18254. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50733-2, EIF2AK3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68482-2-PBS. Detection antibody: 68482-4-PBS. Standard:Ag18254. Range: 0.391-100 ng/mL.