

For Research Use Only

# EIF2AK3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68482-4-PBS



## Basic Information

Catalog Number:

68482-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18254

GenBank Accession Number:

BC126354

GeneID (NCBI):

9451

UNIPROT ID:

Q9NZJ5

Full Name:

eukaryotic translation initiation

factor 2-alpha kinase 3

Calculated MW:

1116 aa, 125 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E12H12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50733-2, EIF2AK3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68482-2-PBS. Detection antibody: 68482-4-PBS. Standard: Ag18254. Range: 0.391-100 ng/mL.