

For Research Use Only

TIA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68486-3-PBS



Basic Information

Catalog Number: 68486-3-PBS	GenBank Accession Number: BC015944	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 7072	CloneNo.: 1G11A4
Source: Mouse	UNIPROT ID: P31483	
Isotype: IgG1	Full Name: TIA1 cytotoxic granule-associated RNA binding protein	
Immunogen Catalog Number: AG2778	Calculated MW: 214 aa, 24 kDa, 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51061-2, TIA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68486-3-PBS. Detection antibody: 68486-1-PBS. Standard:Ag2778. Range: 0.391-50 ng/mL