

For Research Use Only

# FLAD1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68491-4-PBS



## Basic Information

Catalog Number:

68491-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6845

GenBank Accession Number:

BC011378

GeneID (NCBI):

80308

UNIPROT ID:

Q8NFF5

Full Name:

FAD1 flavin adenine dinucleotide  
synthetase homolog (S. cerevisiae)

Calculated MW:

446 aa, 49 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B5C8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50710-2, FLAD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68491-4-PBS. Detection antibody: 68491-3-PBS. Standard: Ag6845. Range: 0.098-100 ng/mL.