

For Research Use Only

# CIAPIN1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68492-4-PBS



## Basic Information

Catalog Number:

68492-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33193

GenBank Accession Number:

BC024196

GeneID (NCBI):

57019

UNIPROT ID:

Q6FI81

Full Name:

cytokine induced apoptosis inhibitor

1

Calculated MW:

312 aa, 34 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2G6B8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50561-2, CIAPIN1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68492-2-PBS. Detection antibody: 68492-4-PBS. Standard: Ag33193. Range: 0.391-100 ng/mL.