

For Research Use Only

UBXN1 Monoclonal antibody, PBS Only



Catalog Number: 68499-1-PBS

Basic Information

Catalog Number: 68499-1-PBS	GenBank Accession Number: BC000902	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 51035	CloneNo.: 3A6C2
Source: Mouse	UNIPROT ID: Q04323	
Isotype: IgG1	Full Name: UBX domain protein 1	
Immunogen Catalog Number: AG33606	Calculated MW: 33 kDa	
	Observed MW: 40 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse, rat

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68499-1-Ig (UBXL1 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 68499-1-PBS in a different storage buffer formulation.