

For Research Use Only

TEX264 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68507-3-PBS



Basic Information

Catalog Number: 68507-3-PBS	GenBank Accession Number: BC008742	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 51368	CloneNo.: 1D8H8
Source: Mouse	UNIPROT ID: Q9Y6I9	
Isotype: IgG1	Full Name: testis expressed 264	
Immunogen Catalog Number: AG33710	Calculated MW: 313 aa, 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50870-2, TEX264 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68507-2-PBS. Detection antibody: 68507-3-PBS. Standard:Ag33710. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50870-3, TEX264 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68507-4-PBS. Detection antibody: 68507-3-PBS. Standard:Ag33710. Range: 3.125-100 ng/mL.