

For Research Use Only

LTF Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68509-4-PBS



Basic Information

Catalog Number: 68509-4-PBS	GenBank Accession Number: BC015822	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 4057	CloneNo.: 4F2E6
Source: Mouse	UNIPROT ID: P02788	
Isotype: IgG1	Full Name: lactotransferrin	
Immunogen Catalog Number: AG33645	Calculated MW: 80 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50548-5, LTF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68509-2-PBS. Detection antibody: 68509-4-PBS. Standard: Ag33645. Range: 1.563-100 ng/mL.