

For Research Use Only

# GATM Monoclonal antibody

Catalog Number: 68513-1-Ig



## Basic Information

**Catalog Number:**

68513-1-Ig

**Size:**

1000 ug/ml

**Source:**

Mouse

**Isotype:**

IgG1

**Immunogen Catalog Number:**

AG3754

**GenBank Accession Number:**

BC004141

**GeneID (NCBI):**

2628

**UNIPROT ID:**

P50440

**Full Name:**

glycine amidinotransferase (L-arginine:glycine amidinotransferase)

**Calculated MW:**

423 aa, 48 kDa

**Observed MW:**

48 kDa

**Purification Method:**

Protein G purification

**CloneNo.:**

1D12G5

**Recommended Dilutions:**

WB 1:5000-1:50000

## Applications

**Tested Applications:**

WB, ELISA

**Species Specificity:**

human, mouse, rat, pig

**Positive Controls:**

WB : pig kidney tissue, rat kidney tissue, mouse kidney tissue, pig pancreas tissue, rat pancreas tissue, mouse pancreas tissue

## Background Information

GATM, also known as AGAT, produces ornithine and guanidinoacetate (GAA) from arginine and glycine, and is a key enzyme for creatine synthesis. GATM deficiency in early infancy causes neurodevelopmental delay. GATM is predominantly expressed in kidney and pancreas in adults. Decrease of GATM proteins was observed in proximal tubule damage cases, which links GATM as a general marker for proximal tubule damage. Alternative splicing of GATM generates three isoforms with molecular weights between 44 kDa to 48 kDa.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



pig kidney tissue were subjected to SDS PAGE followed by western blot with 68513-1-Ig (GATM antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.



Various lysates were subjected to SDS PAGE followed by western blot with 68513-1-Ig (GATM antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours.